Ref. PDI 100606

IgA TURBIDIMETRICO

Conf. 40 ML

IVD

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Conf. 40 ML

IVD

This kit allows the quantitative determination of IgA in human serum by turbidimetric immunoassay.

INTRODUCTION

Immunoglobulins (Ig), also called antibodies, are a group of glycoproteins present in the blood and tissue fluids. Are produced by B lymphocytes in their stage of terminal differentiation of plasma cells or, when the immune system of the subject comes into contact with foreign substances (antigens), belonging for example to bacteria and viruses. The immunoglobulins produced following such stimulation are specific to that antigen and are able to bind with it, facilitating the killing and elimination. Are protein molecules consist of two heavy chains (H, from the English "heavy"), of about 400 amino acids, and two light (L, from the English "light"), of about 200 amino acids held together by chemical bonds. At one end of the antibody binds the antigen. The remaining end performs other functions, such as binding to macrophages (cells responsible for the elimination of foreign substances) and the activation of the complement system.

Classes and functions of immunoglobulins. Are detectable five classes of immunoglobulins: IgA, IgD, IgE, IgG and IgM, each of which possesses a specific class of H chains, which are respectively alpha, delta, epsilon, gamma and mu. The IgA antibodies are more present in secretions: saliva, milk, tears, and respiratory secretions, digestive and genitourinary. Can exist in both single units (monomers), or aggregated in pairs (dimers). They are a real defensive system in the mucous membranes, representing the first specific barrier that opposes the penetration of antigenic material in the body. The determination of IgA antibodies is important to differentiate immune deficiencies and myelomas. It also plays a role in acute and chronic infections as a first line of defense.

Increased levels can be found in acute infectious hepatitis, chronic aggressive hepatitis, alcoholic cirrhosis, active chronic infection, rheumatoid arthritis, mixed connective tissue disease.

PRINCIPLE

The reagent determines the concentration of IgA present in the serum sample against a standard curve to known value. The addition of serum containing IgA to a solution that contains the corresponding antiserum, determines a turbidity that allows, by the OD value obtained, to trace the concentration of antigen present in the sample.

REAGENTS

For in vitro diagnostic use.

The reagent stored at 2-8°C, is stable until the expiration data printed on the label. Do not freeze the kit components. The liquid reagent is ready to use.

Kit Components:		
Reagent A	2x20 ml	Ref.
Human serum Anti IgA		PDI 1006061
Polyclonal goat anti-human IgA stabilized in saline solution with integrated accelerator.		
Sodium azide (0.95 g / L).		
System Components supplied separately		
Low and High Control	1x1 ml	
Pooled human sera positive for Antithrombin III, IgA, IgG, IgM, and transferrin. The value is stated on the label.		
High Calibrator	1x1 ml	
or		
Calibrators Set	5x1 ml	
Pooled human serum calibrated to the Antithrombin III, IgA, IgG, IgM, and transferrin. The value is stated on the label.		

SAMPLE

Use fresh serum. If the test can not be performed on the same day, the serum may be stored at 2-8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

PROCEDURE

Manual Procedure

Sample/Control: diluite 1:10 with saline solution 9 g/L
Reference curve: create a reference curve by diluting the high level of the standard 1:10, 1:20, 1:40, 1:80 and 1:160 with a saline 9 g/L. Use a saline 9 g /L as the zero point.

Note: If in possession of a set of calibrators to 5 levels, use these standards as they are to create the curve.

Test: mixing 16 µl of the diluted samples, standards and control (s) with 1000 l of Reagent A and incubate for 5 minutes at room temperature. Read the optical density (OD) of the samples, standards and control (s) at 340 nm. Draw a standard curve and read the concentration of samples and controls, taking into account in the final calculation of the dilution performed on the sample/control (1:10).

AUTOMATION

The application of automatic analyzers such as Hitachi and PKL 300 PPC are available on request.

REFERENCE VALUES

Man: 83 - 406 mg/dL (IFCC)

Woman: 70 - 374 mg/dL

These values are indicative. Each laboratory should establish its own reference values.

PERFORMANCE

The performance for the monoreagent Turbidimetric IgA were measured on a clinical chemistry analyzer ³

(Selectra 2).

Range of measurement: 0-654 mg / dL

Limit of sensitivity: 20 mg / dL

Hook effect: > 13080 mg/dL

Sensitivity: 0,00045 unit ABS / concentration unit

Precision:

[%CV]

	Low	Medium	High
Intra-Run	1.72	2.28	1.59
Inter-Run	3.34	2.97	1.96

Accuracy:

[mg/dL]

Control	attached	measured
BEHRING	160 (136 – 184)	146
BIORAD 1	126 (107 – 145)	128

Specificity:

monospecific

Interferences: No interference by: hemoglobin (1000 mg/dL), Na-citrate (1000 mg/dL), heparin (50 mg/dL), triglycerides (2500 mg/dL), EDTA (5mg/dl). Turbidity (>2.5%) and bilirubin (> 5mg/dl) interfere with the test.

Limits:none

Compare another similar commercially available reagent:
y = 1.3167x - 67,272 / r = 0,9913

Stability at 2 - 8 °C: at least 3 years after production

PRECAUTIONS and WARNINGS

- 1.Only for in vitro diagnostic use.
- 2.Sodium azide is able to form lead or copper azides in laboratory plumbing which may explode. Rinse thoroughly with water after disposing of waste fluids containing sodium azide.
- 3.Each donor unit used in the preparation of standards and controls were negative for the presence of antibodies to HIV1 and HIV2, as well as for hepatitis B surface antigen and hepatitis C antibodies, using a method approved by the FDA.

WASTE MANAGEMENT

Please refer to official regulations.

CLASSIFICATION EDMA

Name IgA ; n°: 12.01.01.01.00

BIBLIOGRAPHY

- 1.Dati, F. et al., Lab. Med. 13, 87 (1989)

Ver. 001 del 15.10.2012