Ref. PDI103410

CK-NAC UV LIQUIDA

Conf. 125 ml

IVD

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Nome CK-NAC n:11.01.01.13.00.

BIBLIOGRAFIA

Methods of Enzymatic Analysis, (1983), 3th edition, V III Edited by H. U. , Bergmeyer, Verlag Chemie.

Ref. PDI103410

CK-NAC UV LIQUID

Conf. 125 ml

IVD

INTRODUCTION

CPK is an intracellular enzyme found mainly in skeletal muscle, myocardium and brain tissue. Catalyzes the phosphorylation of creatine, thanks to the intervention of ATP, with formation of phosphocreatine and ADP, reaction of fundamental importance in muscle metabolism. Have been identified in three different isoenzymes diagnostic significant: CPK-BB of cerebral origin, CPK-MB of cardiac origin and CPK-MM of musculoskeletal origin. Exercise, deep intramuscular injections, delirium tremens, and surgery may raise the levels of CPK. Drugs that may produce an increase of CPK are anticoagulants, morphine, alcohol, high doses of salicylates, amphotericin B, clofibrati and some anesthetics. It is advisable to pick fasted for 8 hours.

PRINCIPLE

IFCC method. The creatine kinase catalyzes the reaction between creatine phosphate and ADP, with formation of creatine and ATP, which, in the presence of glucose and hexokinase (HK), is transformed into ADP and glucose-6-phosphate. The latter, in the presence of glucose 6-phosphate dehydrogenase (G-6-PDH), reacts with b-NADP + to form 6-phosphogluconate and b NADPH. By measuring the change in absorbance due to the conversion of b- NADP+ in b-NADPH, in the unit of time at 340 nm is calculated CK activity in the sample. The presence of N-acetylcysteine (NAC) in the reaction mixture allows for optimal reactivation of the enzyme CK.

METHOD

IFCC.

REAGENTS

For in vitro diagnostic use. The reagents, stored at 2-8°C, are stable until the expiration date printed on the label.

Kit Components:

A) Reagent A: buffer	1x100 ml	Ref. PDI 1034101
B) Reagent B: substrate	1x25 ml	Ref. PDI 1034102

Monoreagent PREPARATION

(if you work with monoreagent procedure)

Mix one part of Reagent 2 (Liquid Substrate) with 4 parts of Reagent 1 (buffer). This solution remains stable 2 weeks, when stored in the refrigerator at 2-8 °C or 2 days when stored at room temperature (<25 °C). Protect from direct sunlight.

SAMPLE

Serum or plasma (heparin or EDTA). Hemolysis up to 200 mg Hb/dl does not interfere in the determination.

Concentration of each substance in the final reaction volume:

Imidazole acetate buffer	100 mmol/l pH 6.8
creatine phosphate	50 mmol/l
D-glucose	20 mmol/l
Mg acetate	10 mmol/l
NAC	20 mmol/l
EDTA	1.6 mmol/l
ADP	3 mmol/l
AMP	5 mmol/l
NADP	2 mmol/l
hexokinase	>5000 U/l
diadenosine pentaphosphate	20 µmol/L
G-6-PDH	>7000 U/L

Preservatives and stabilizers.

Warning: Contains sodium azide (0.1%) as preservative.

OPERATING CONDITIONS

Wavelength:	340 nm
Optical path:	1 cm
Reaction Temperature:	37°C
Reaction:	kinetic

Allow reagents to working temperature before use.

ANALYTICAL PROCEDURE

Monoreagent Procedure

Pipette in a thermostated cuvette	
Working solution	200 µL
Allow reagents to working temperature before use (2-3 minutes at 37 °C)	
Sample	8 µL
Mix and take the first reading after 100 seconds. Repeat the reading of more times in a time interval of 200 seconds. Calculate the average D A /min found.	

ANALYTICAL PROCEDURE

Bireagent Procedure

Mix and take the first reading after 100 seconds. Repeat the reading of more times in a time interval of 200 seconds. Calculate the average D A / min found.	
Pipette in a thermostated cuvette	
Reagent 1 buffer	200 µL
Sample	10 µL
Mix, incubate for about 1 minutes and then add:	
Reagent 2 Substrate	50 µL
Mix and take the first reading after 120 seconds. Repeat the reading of more times in a time interval of 200 seconds. Calculate the average D A / min. found.	

CALCULATION

Multiply the DA/min average found by the factor 4127 valid for the operating conditions.

STANDARD VALUE

	37°C
Man:	24-170 U/l
Woman:	24-140 U/l

Activities for AE/min. above 0. 100 (365nm) and to 0. 180 (340-334nm), repeat the determination on the sample diluted 1:10 with saline and then multiply the result by 10.

CALIBRATORS AND/OR CONTROLS

For internal quality control using PKL Normal Control Serum Human or PKL Pathological Control Serum Human.

It is advisable to analyze the two levels of checks at least once a day. Values obtained should fall within the specified range.

If the values should be outside this range and if any error in question was excluded from retesting, you should apply the following points:

1. Check the programming of the instrument and the light source.
2. Check cleanliness of all equipment in use.
3. Check the water contamination (eg. Bacterial) may contribute to inaccurate results.
4. Check the reaction temperature.
5. Check the expiry date of the kit and its contents.

ANCILLARY EQUIPMENT

Solution of NaCl 9 g/l. Various volume pipettes, spectrophotometer or photometer Laboratory.

PRECAUTIONS

The reagent contains sodium azide as preservative (0.1%). Do not swallow! Avoid contact with skin and mucous membranes. In case of accidental spillage dilute with water. **These products are for professional use only and should be used by qualified personnel.**

SENSIBILITY

The kit determines a minimum concentration of 6 U/L of CK-NAC in serum.

INTERFERENCES

Bilirubin concentrations up to 12 mg/dl do not cause interferences.

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