Ref. PDI 125060

CREATININA ENZIMATICA LIQUIDA

Conf. 100 ml

IVD

INDICAZIONE

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Nome Creatinina n. 11.02.01.07.00.

BIBLIOGRAFIA

Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982). National Kidney Foundation K/DOQI. Clinical Practice Guide-lines for chronic kidney disease: evaluation, classification, and stratification. Am J Kidney Dis 2002; 39:S1-S200

Badiou S, Dupuy AM, Descomps B, Cristolead, JP. Comparison between the enzymatic vitros assay for creatinine determination and three other methods adapted on the Olympus analyzer, Journal of Clinical Laboratory Analysis 2003;17, 235 - 240. Hayes AW. Principles and Methods of Toxicology, Taylor & Francis, 1028 (2001)

Cristenson RH, Johnson LJ, Gregory, LC. Appleton and Lange's Outline Review Clinical Chemistry, McGraw-Hill Professional, 118 (2001). Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatinine Assays in Plasma and Serum and Early Morning Urine. Clin Lab 2000; 46;53-55.

Ref. PDI 125060

CREATININE ENZYMATIC LIQUID

Conf. 100 ml


IVD

INDICATION

The creatinine is derived from creatine and from creatine phosphate in muscle tissues and can be defined nitrogenous waste product. Creatinine is not reused and is excreted from the body in the urine through the kidneys. It is produced and excreted at a constant speed, speed that is proportional to the muscle mass body. Being excreted through the kidneys, the measurement of creatinine is used almost exclusively in the monitoring of renal function. It is considered as the most useful endogenous marker in the diagnosis and treatment of renal diseases. The measurement of creatinine instead of urea for controlling the operation renal presents some advantages. The plasma creatinine levels tend not to be affected by the ingestion of protein and water, the speed of production of urine and exercise. Having a production rate constant, elevated serum creatinine are indicative of decreased glomerular filtration and suggest kidney dysfunction. Low values of serum creatinine are rare and not significant from a clinical point of view.

PRINCIPLE

For enzymatic colorimetric determination of creatinine in serum, plasma and urine. The method does not require deproteinization.

LEVY AND PREPARATION SAMPLE

Serum or plasma (lithium heparin).
24 hours urine: dilute 1:10 (1 + 9) con saline.

KIT COMPONENT

1. REAGENT 1 (Ref. PDI1250601) 5x15 ml
Creatinase 12 IU-60 IU/ml
TOOS 0.07 mg-0.21 mg/ml
SOD 4-17 IU/ml

2. REAGENT 2 (Ref. PDI1250602) 2x12.5 ml
Creatinase 135 IU ~ 670 IU/ml
4-AA (4-aminoantipirina) 0.3 – 0.9 mg/ml
POD ≥ 20 KU/L

All reagents are ready to use.

WARNINGS AND PRECAUTIONS

Only for in vitro diagnostic use. The product contains active components such as preservatives (sodium azide and other) and detergents. The total concentration of these components is less than the limits given in Directives 67/548/EEC and 1999/45/EC and its amendments and additions to the classification, packaging and labeling of dangerous preparations. It is recommended to not pipette by mouth pipetting. Comply with the rules in force for handling laboratory reagents. The reagents must be used by qualified laboratory technicians in laboratory conditions and appropriate only for the stated purpose.

STABILITY AND PREPARATION OF REAGENTS

The reagents when stored at 2-8 °C, are stable until the expiration date. The stability persists even after the opening of the bottles, if not contaminated during use and if closed immediately after use.

NOTES

A slight yellow coloration of the reactants does not involve any change in the performance of the product.

CALIBRATIONS AND / OR CONTROLS

For calibration use PKL calibrator. For internal quality control using PKL Normal or Pathological Human Control Serum. It is advisable to analyze the two levels of checks at least once a day. Values obtained should fall within the specified range. **If the values should be outside this range and if any error in question was excluded from retesting, you should apply the following points:**

1. Check the instrument settings and the source of light.
2. Check cleanliness of all equipment in use.
3. Check water; any contamination (eg. Bacterial) may contribute to inaccurate results.
4. Check the reaction temperature.
5. Check the expiry date of the kit and its contents.

STANDARD VALUE

Serum:	Men	59 - 104 µmol/l (0.67 - 1.17 mg/dl)
	Women	45 - 84 µmol/l (0.51 - 0.95 mg/dl)
Urine:	Men	3540-24600 µmol/l (40-278 mg/dl)
	Women	2550-20000 µmol/l (29-226 mg/dl)

It is recommended that each laboratory establish its own reference range.

INTERFERENCES

No influence from hemoglobin <1000 mg/dl, triglycerides <1000 mg/dL, bilirubin <40 mg/dl, ascorbic acid <10 mM.

MEASURING RANGE

Serum: 0.1 – 13.5 mg/dl. Urine: 0.1 – 13.5 mg/dl.

For higher concentrations, repeat the determination of a sample diluted 1:10 with saline and then multiply the result by 10.

SENSIBILITY

0.20 mg/dl. The sensibility was calculated on 20 samples of saline and expressed as the "average value of zero + DS 2576."

PERFORMANCE (on automatic Hitachi 911 analyzer)

Serum

Serum Testing	Within-Run Precision			
	0.75 mg/dL (66.3 μM)	1.41 mg/dL (125 μM)	4.11 mg/dL (346 μM)	10.28 mg/dL (908.7 μM)
No. of Data Points	80	80	80	80
Mean mg/dL (μM)	0.74 (65.4)	1.38 (122.3)	4.04 (357.5)	10.28 (908.7)
SD mg/dL (μM)	0.015 (1.3)	0.015(1.37)	0.029(2.54)	0.015 (1.3)
Cv%	2.1	1.1	0.7	0.1
Serum Testing	Total Precision			
	0.75 mg/dL (66.3 μM)	1.41 mg/dL (125 μM)	4.11 mg/dL (346 μM)	10.28 mg/dL (908.7 μM)
No. of Data Points	80	80	80	80
Mean mg/dL (μM)	0.74 (65.4)	1.38 (122.3)	4.04 (357.5)	10.28 (908.7)
SD mg/dL (μM)	0.022(1.9)	0.026(2.29)	0.058(5.11)	0.14(12.4)
Cv%	3.0	1.9	1.4	1.4

Cristenson RH, Johnson LJ, Gregory, LC. Appleton and Lange's Outline Review Clinical Chemistry, McGraw-Hill Professional, 118 (2001). Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatinine Assays in Plasma and Serum and Early Morning Urine. Clin Lab 2000; 46:53-55.

Ver. 001 del 20.05.2010

Urine

Urine Testing	Within-Run Precision		
	Level 1	Level 2	Level 3
No. of Data Points	21	21	21
Mean mg/dL (μM)	29.09 (2572)	87.1 (7711)	196.7 (17407)
SD mg/dL (μM)	0.1 (8.84)	0.27 (23.60)	0.90 (79.71)
Cv%	0.36	0.31	0.46
Urine Testing	Total Precision		
	Level1	Level 2	Level 3
No. of Data Points	20	20	20
Mean mg/dL (μM)	29.86 (2640)	87.7 (7765)	195 (17265)
SD mg/dL (μM)	0.79 (69.8)	0.67 (59.2)	1.19 (105.2)
Cv%	2.64	0.76	0.60

CORRELATION

The correlation between this test (y) with other commercial kit (x) gave the following results:

SERUM: Correlation Coefficient: 0.9981
Slope/Intercept: $y = 0.9467x + 0.0643$

URINE: Correlation Coefficient: 0.9969
Slope/Intercept: $y = 1.005x - 0.2979$

CLASSIFICATION EDMA

Name Creatinine n. 11.02.01.07.00.

BIBLIOGRAPHY

Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982). National Kidney Foundation K/DOQI. Clinical Practice Guide-lines for chronic kidney disease: evaluation, classification, and stratification. Am J Kidney Dis 2002; 39:S1-S200

Badiou S, Dupuy AM, Descomps B, Cristolead, JP. Comparison between the enzymatic vitros assay for creatinine determination and three other methods adapted on the Olympus analyzer, Journal of Clinical Laboratory Analysis 2003:17, 235 - 240. Hayes AW. Principles and Methods of Toxicology, Taylor & Francis, 1028 (2001)

PKL PPC 125 ANALYZER APPLICATION

<p>CUSTOMER</p> <p>PARAMETER</p> <p>QC</p> <p>SCHEDULE</p> <p>REPORT</p> <p>STATISTICS</p> <p>MAINTENANCE</p> <p>RUN MONITOR</p> <p>EXIT</p>	<p>Test Parameter Profile Item-Test sequence Calculation item External parameter Reflex test</p> <p>Full name CREE Item CREE Text Code XXX</p> <p>Basic parameter Reference range Calibration</p> <p>Method Two points Primary filter 546 Secondary filter NO</p> <p>Reagent blank Reagent blank Decimal 2 Unit Mg/dl</p> <p>R1 Volume 225 R1 Position 7 Incubation time (s) 210</p> <p>R2 Volume 75 R2 Position 8 Incubation time (s) 40</p> <p>Sample volume 7 Read time (s) 144 Max speed set middle</p> <p>Reagent suppliers paramedical Reagent barcode</p> <p>Lot number Expiry Date Lot number</p> <p> <input type="checkbox"/> Reagent blank 0.0000 Lower 0.0000 High 0.3450 Co-relation Y= 1.0 X+ 0.0 Dilution ratio Linear range 1000.00 </p> <p> <input type="button" value="Add"/> <input type="button" value="Delete"/> <input type="button" value="Save"/> <input type="button" value="Printview"/> <input type="button" value="Print"/> </p>
--	--

XXX = Entered by operator