Ref. PDI 125220

α -1 Glicoproteina acida

Conf. 25 ml

IVD

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Nome: α -1 Glicoproteina acida ; n°: 12.01.90.01.00

BIBLIOGRAFIA

1. Schmid, K. in FW putman, Editor, The plasma Proteins, Vol 1, second edition, Academic Press, new York, 1975, pp184-228
2. Jonhson, A.M. et al., J. Clin. Invest., **48**(1969) 2293
3. Dati, F. et al., Lab. Med. **13** (1989) 87

Ref. PDI 125220

α -1 Acid Glycoprotein

Conf. 25 ml

IVD

INTRODUCTION

As an early acute phase reactant, AGP is especially useful in monitoring tumor recurrence. Levels are also helpful in differentiating acute phase responses (elevated levels) from estrogen effects (normal or depressed levels). In addition it is an excellent protein to assay along with Haptoglobin in assessing in vivo hemolysis. An elevated AGP level but normal Haptoglobin suggests an acute phase response with mild to moderate in vivo hemolysis.

PRINCIPLE

Measurement of antigen-antibody reaction by the end-point method.

REAGENTS

For in vitro diagnostic use.

The reagents stored at 2-8°C, are stable until the expiration data printed on the label. Do not frozen the kit components.

Kit Components:		
Reagent A	2x12.5 ml	Ref. PDI1252201
Polyclonal goat anti-human AGP stabilised in saline supplemented with accelerator. Sodium azide (0.95 g/L).		
System Components supplied separately		
Standard High Proteins	1x1 ml	
Pooled human sera. The value is stated on the label.		

SAMPLE

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2 – 8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

MANUAL PROCEDURE

Sample/Control: Dilute 1:10 with saline 9g/L.

Calibration curve: generate a calibration curve to 5 points. Use saline solution as zero and the standard as such as the end point (S5). Proceed then with a 1:1 dilution of the standard to determine the other 3 standard (S2, S3, S4) at a known concentration.

STANDARD VALUES

Man: 50 - 130 mg/dL (IFCC)

Woman: 40 - 120 mg/dL

This range is given for orientation only. Each laboratory should establish its own reference values.

PERFORMANCES

The performance characteristics for the AGP were measured on a clinical chemistry analyzer (Selectra 2).

Measuring Range: 0 - 299 mg/dL

Detection Limit: 9 mg/dL

Hookeffect: > 400 mg/dL

Sensitivity: 0.00158 ABS units/concentration unit

Precision:		Low	Medium	High
[%CV]	Intra-Run	1.29	1.30	1.30
	Inter-Run	2.23	2.08	2.17

Accuracy:	Control	Assigned	Measured
[mg/dL]	BEHRING	71.7 (60.9 – 82.5)	63.4
	CLINIQA III	128 (109.2 – 147.2)	128.5

Specificity: Monospecific

Interferences: No interference for: Na-citrate (1000 mg/dL), Heparin (50 mg/dL), Triglyceride (2500 mg/dL) and EDTA (5 mg/dL). Haemoglobin (250 mg/dL), Turbidity (> 2.5 %) and Bilirubin (> 10 mg/dL) interfere with the test.

Limitations: None

Comparison with Bireagents kit:

$y = 0.9141x - 0.8383 / r = 0.9674$

Stability at 4° C: at least 3 years after production

PRECAUTIONS and WARNINGS

1. In vitro diagnostic use only.
2. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.
3. Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

CLASSIFICATION EDMA

Name: α -1 Acid Glycoprotein ; n°: 12.01.90.01.00

BIBLIOGRAPHY

4. Schmid, K. in FW putman, Editor, The plasma Proteins, Vol 1, second edition, Academic Press, new York, 1975, pp184-228
5. Jonhson, A.M. et al., J. Clin. Invest., 48 (1969) 2293
6. Dati, F. et al., Lab. Med. 13 (1989) 87

Ver. 002 del 07.02.2010

CUSTOMER

PARAMETER

QC.

SCHEDULE

REPORT

STATISTICS

MAINTENANCE

RUN MONITOR

EXIT


Test Parameter
Profile
Item-Test sequence
Calculation item
External parameter
Reflex test

Full name

Item

Test Code

Basic parameter
Reference range
Calibration

K =

Number of standards

Calib Curve

Calibration cup

Cali. Hist.

1

Standard position

Standard value

Absorbance

Lot number

Expiry Date