Ref. PDI125230

BILIRUBINA DIRETTA

Conf. 110 ml

IVD

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Nome Bilirubina Diretta n. 11.01.02.03.00.

BIBLIOGRAFIA

Charles T.W. Lam, Elkhart, Test composition, device and method for detecting Bilirubin. (1975).

Ferrero C, Modenese A, Bonini PA, Ceriotti F, Carobene A, Franzin M, G. Guder W, Hofmann W, PUUKKA M, Laitinem H, Makarovsy V and Gaidukova N. Multicentre Evaluation of KONE Optima Analysis System. Clin Chem Lab Med 1998; 36(7): 475-484.

Martin H. Kroll, Ronald J. Elin. Interference with Clinical Laboratory Analyses. Clin. Chem. 40/11(1994).

Philippe Marquis. Comparison of analytical methods in clinical chemistry. Service de biochimie Centre Hospitalier, Metez-France.

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IVD

INDICATION

Bilirubin is a breakdown product of hemoglobin. Free, unconjugated, is nearly apolar and practically insoluble in water, bound to albumin forms a soluble complex that can transport blood from the spleen to the liver. In the liver, bilirubin is conjugated with glucuronic acid and the resulting product is soluble in water, is excreted through the bile ducts. Hyperbilirubinemia may be caused by increased bilirubin due to hemolysis (pre-hepatic jaundice), damage to liver parenchyma (intra-hepatic jaundice) or occlusion of the bile ducts (post-hepatic jaundice). A chronic form of congenital hyperbilirubinemia (predominantly unconjugated form), called Gilbert syndrome is very common in the population. Elevated levels of total bilirubin were observed in 60-70% of infants after increased degradation of erythrocytes after birth and for the delayed operation of the enzymes used for degradation of bilirubin. Common methods of determination of total bilirubin bilirubin that is dosed on the air. The determination of direct bilirubin measure predominantly conjugated bilirubin, soluble in water. The unconjugated bilirubin can therefore be estimated from the difference between the total and direct bilirubin.

ASSAY PRINCIPLE

The test is based on a method modified Jendrassik-Grof. The direct bilirubin, in acid solution, reacts with diazotized sulfanilic acid to form a colored azocompound (azobilirubin), whose intensity is directly proportional to the concentration of direct bilirubin in the sample.

COMPOSITION

A) REAGENT A: 5x20 ml (Ref. PDI1252301)
Sulfanilic Acid 30 mmol/l
Hydrochloric Acid 170 mmol/l

B) REAGENT B: 1x10 ml (Ref. PDI1252302)
Sodium Nitrite 29 mmol/l

The reagents are ready for use.

CONSERVATION AND STABILITY

Store at room temperature (18-25 °C).

Do not freeze the reagents!

After opening:

- REAGENT A and REAGENT B are stable until the expiration date printed on the label, if you avoid contamination and evaporation.
- The above conditions are valid if the bottles are opened only for taking the time to take the desired amount, immediately closed with its cap and stored at the temperature indicated.

ANCILLARY EQUIPMENT

- Automatic pipettes
- Photometer
- Sample cuvette (optical path = 1 cm)
- Thermostatic cell
- NaCl solution (9 g/l)
- Calibrator Bilirubine or Multicalibrator

SAMPLE

Fresh serum not hemolyzed or plasma. Serum should be analyzed within 2 hours. The turbidity of the sample due to the presence of macromolecules may interfere. In this case, we suggest the centrifugation or filtration through membrane filters of 0.2 µ. Avoid exposure to light.

STANDARD VALUE

Adult: up to 0.25 mg/dl (4.3 µmol/L).

Each laboratory should establish its own normal values according to the population in which it operates.

SENSIBILITY

The sensibility of the method (Average + 3 DS), is of 0.03 mg/dl (0.51 µmol/L).

MEASURING

The test was developed for determination of direct bilirubin in a measuring range of 0.03 - 20 mg/dl (0.51 - 342 µmol/L). If these values are exceeded, samples should be diluted 1 +9 with a NaCl solution (9 g/l) and the result multiplied by 10.

CORRELATION

The comparison of this method(y) with a similar method(x) gave the following correlation:

$$y = 1.0232x - 0.7022 \quad r = 0.9984 \quad N = 50$$

NOTE

The reaction is inhibited by reducing agents, samples of patients with massive ascorbic acid may cause falsely low.

CALIBRATOR AND/OR CONTROL

It is recommended in the execution of the method, type of control sera with values of direct bilirubin are known to check compliance of data and planned results and thus validate the data. It is advisable to analyze two levels of checks at least once a day. The values obtained must be in the range specified. If values were to be outside this range and if a possible error in question was excluded from the repetition of the test, you should apply the following points:

1. Check instrument settings and the source of light.
2. Check cleanliness of all equipment in use.
3. Check water, possible contamination (eg. Bacterial) may contribute to inaccurate results.

4. Check the temperature of reaction.
5. Check the expiry date of kit and its contents.

WASTE MANAGEMENT

Please refer to local legal requirements.

PRECAUTIONS

The reagent contains sodium azide as preservative (0.1 %).

Do not swallow! Avoid contact with skin and mucous membranes. In case of spillage, dilute with water.

These products are only for professional use and should be used by qualified personnel.

SPECIFICITY

The reagent only determines Direct Bilirubin in samples of human origin.

CLASSIFICATION EDMA

Name Direct Bilirubin n. 11.01.02.03.00

Ver. 002 del 07.07.2010

BIBLIOGRAPHY

Charles T.W. Lam, Elkhart, Test composition, device and method for detecting Bilirubin. (1975).

Ferrero C, Modenese A, Bonini PA, Ceriotti F, Carobene A, Franzin M, G. Guder W, Hofmann W, PUUKKA M, Laitinen H, Makarovsy V and Gaidukova N. Multicentre Evaluation of KONE Optima Analysis System. Clin Chem Lab Med 1998; 36(7): 475-484.

Martin H. Kroll, Ronald J. Elin. Interference with Clinical Laboratory Analyses. Clin. Chem. 40/11(1994).

Philippe Marquis. Comparison of analytical methods in clinical chemistry. Service de biochimie Centre Hospitalier, Metz-France.

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CUSTOMER
PARAMETER
QC
SCHEDULE
REPORT
STATISTICS
MAINTENANCE
RUN MONITOR
EXIT

Test Parameter
Profile
Item-Test sequence
Calculation item
External parameter
Reflex test

Full name:
Item:
Test Code:

Basic parameter
Reference range
Calibration

K =

Number of standards:
Calib Curve:
Calibration cup:
Cali. Hist.

1
Standard position:
Standard value:
Absorbance:

Lot number:
Expiry Date:

Add
Delete
Save
Printview
Print

XXX = Entered by the operator