Ref. PDI 125310

**RAME**

Conf. 60 ml

**IVD**

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Nome Rame n.: 11.02.01.06.00.

## BIBLIOGRAFIA

Joseph C. Robinson. A simple method for determining serum copper (1949).

Siedel, Dr. Joachim , Schellong, Lieselotte, Staepels, Johnny, Hermann, Dr. Uwe, Town, Dr. Michael-harold Process and reagent for the determination of iron. (1988).

Ferrero C, Modenese A, Bonini PA, Ceriotti F, Carobene A, Franzin M, G. Guder W, Hofmann W, Puukka M, Laitinen H, Makarovsy V and Gaidukova N. Multicentre Evaluation of KONE Optima Analysis System. Clin Chem Lab Med 1998; 36(7): 475-484.

Bablok W, Passing H. Application of statistical procedures in analytical instrument testing. J Autom Chem 1985; 7:74-9.

Martin H. Kroll, Ronald J. Elin. Interference with Clinical Laboratory Analyses. Clin. Chem. 40/11 (1994).

Philippe Marquis. Comparison of analytical methods in clinical chemistry. Service de biochimie Centre Hospitalier, Metz- France.

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IVD

### INDICATION

Copper is a co-factor of different enzymes such as cytochrome oxidase, tyrosinase, uricase. It is involved in iron metabolism, promoting its intestinal uptake, and its mobility from storage tissue. A decrease level of copper is correlated to a decreased level of proteins in the serum, therefore in case of inadequate nutrition or not correct uptake (celiac disease, sprue), loss of protein through faces or urine (nephritic syndrome), or in Wilson disease. An increase of copper level occurs in pregnancy, acute or chronic infections, surgery, myocardium infarction, iperthyroidism, and haematological diseases.

### ASSAY PRINCIPLE

Copper (Cu<sup>++</sup>) reacts with the chromogen Di-Br-PAESA at room temperature yielding a coloured complex which intensity is proportional to the Copper concentration present in the sample. The method does not require serum deproteinisation either sample blank.

### REAGENTS

For in vitro diagnostics use.

Kit components:

| Reagent A                                                                                     | 4x12.5 ml | Ref. PDI 1253101 |
|-----------------------------------------------------------------------------------------------|-----------|------------------|
| Ready for use. Stored at 2-8°C remains stable until the expiration date printed on the label. |           |                  |
| Acetate buffer pH 4.7                                                                         |           | 0.160 M          |
| Surfactants                                                                                   |           | 1.5%             |
| Preservatives and stabilizers                                                                 |           |                  |
| Reagent B                                                                                     | 1x10 ml   | Ref. PDI 1253102 |
| 3,5-dibromo-2-pyridylazo<br>N-ethyl-N-3-sulfopropyl-aniline                                   |           | 0,1 mM           |
| Preservatives and stabilizers                                                                 |           |                  |

### SAMPLES

Use fresh serum or plasma without hemolysis.

### ANALYTICAL PROCEDURE

#### ATTENTION!!!

**Reagent A may present precipitates when stored at low temperatures. In this case heated to 37°C before use for the disappearance of these precipitates.**

### STANDARD VALUES (in serum)

Man: 70 µg/dl - 140 µg/dl (11.0-22.0 µmol/L)

Women: 80 µg/dl - 155 µg/dl (12.6-24.4 µmol/L)

It is recommended that each laboratory establish its own expected range.

### LINEARITY

The reaction is linear up to 400 µg/dl. If you get higher concentrations, repeat the test serum diluted 1 +1 with saline, then multiplying the result by 2.

### CALIBRATOR AND/OR CONTROL

For calibration using the PKL calibrator. For the internal quality control using PKL control sera of normal human or the PKL human control serum Pathological.

It is advisable to analyze two levels of checks at least once a day. The values obtained must be in the range specified. If values were to be outside this range and if a possible error in question was excluded from the repetition of the test, you should apply the following points:

1. Check instrument settings and the source of light.
2. Check cleanliness of all equipment in use.
3. Check water, possible contamination (eg. Bacterial) may contribute to inaccurate results.
4. Check the temperature of reaction
5. Check the expiry date of kit and its content.

### INTERFERENCES

Do not use hemolyzed serum samples. It should use disposable plastic tubes to avoid any possible contamination and interference.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### ANCILLARY EQUIPMENT

NaCl solution 9 g/l. Various volume pipettes, spectrophotometer or photometer Laboratory.

### PRECAUTIONS

The reagent contains sodium azide as preservative (0.09 %). Do not swallow! Avoid contact with skin and mucous membranes. In case of spillage, dilute with water.

**This products are for professional use only and should be used by qualified personnel.**

### SPECIFICITY

The reagent determines only the copper in human samples.

### SENSIBILITY

The minimum detectable activity is 2 µg/dl with a Hitachi 911 Analyzer.

### ACCURACY

The comparison of this test copper (y) and the atomic absorption method (x) gave the following results:

N = 60, r = 0.9612, y = 0.9771 x + 2.5769

### PRECISION

| Serum Control | Lot    | Value Obtained | Average | Range    | Recovery % |
|---------------|--------|----------------|---------|----------|------------|
| PKL Normal    | 228 UN | 98             | 100     | 80 - 119 | 102        |

|                     |        |     |     |           |      |
|---------------------|--------|-----|-----|-----------|------|
| PKL<br>Pathological | 141 UE | 171 | 168 | 143 - 193 | 98.2 |
|---------------------|--------|-----|-----|-----------|------|

### REPEATABILITY AND REPRODUCIBILITY

The repeatability is expressed as a coefficient of variation in series.  
The reproducibility is expressed as a coefficient of variation between sets.

| CV (%)       | Level I<br>33 µg/dl | Level II<br>88 µg/dl | Level III<br>219 µg/dl |
|--------------|---------------------|----------------------|------------------------|
| Intra series | 4.41                | 3.09                 | 2.24                   |
| Inter series | 5.03                | 3.82                 | 2.03                   |

Ver. 001 del 10.06.2010

### CLASSIFICATION EDMA

Name Copper n. 11.02.01.06.00.

### BIBLIOGRAPHY

Joseph C. Robinson. A simple method for determining serum copper (1949).

Siedel, Dr. Joachim , Schellong, Lieselotte, Staepels, Johnny, Hermann, Dr. Uwe, Town, Dr. Michael-harold Process and reagent for the determination of iron. (1988).

Ferrero C, Modenese A, Bonini PA, Ceriotti F, Carobene A, Franzin M, G. Guder W, Hofmann W, Puukka M, Laitinen H, Makarovsy V and Gaidukova N. Multicentre Evaluation of KONE Optima Analysis System. Clin Chem Lab Med 1998; 36(7): 475-484.

Bablok W, Passing H. Application of statistical procedures in analytical instrument testing. J Autom Chem 1985; 7:74-9.

Martin H. Kroll, Ronald J. Elin. Interference with Clinical Laboratory Analyses. Clin. Chem. 40/11 (1994).

Philippe Marquis. Comparison of analytical methods in clinical chemistry. Service de biochimie Centre Hospitalier, Metez- France.

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
Basic parameter
Reference range
Calibration

K =

Number of standards
No
Calib Curve
26.0
Calibration cup
Cup
Cali. Hist.

1

Standard position

Standard value

Absorbance

26.0

26.0

Lot number

Expiry Date

Add
Delete
Save
Printview
Print