# MacCONKEY MUG AGAR

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pH = 7.1 ± 0.2

## PREPARATION

Suspend 46.3 gr of powder in 1 litre of distilled water. Heat until completely dissolved. Sterilize in the autoclave at 121°C for 15 minutes.

## USE

MacCONKEY MUG AGAR is a selective medium recommended for the isolation of enterobacteria, for the differentiation of coliform organisms from lactose non fermenting pathogenic enterobacteria and for the rapid identification of *Escherichia coli*. The colorless MUG (4-methylumbelliferil-β-D-glucuronide) is hydrolyzed by the enzyme β-glucuronidase present in the *E. coli* while releasing the compound 4-methylumbelliferone, highly fluorescent if observed underneath a UV lamp with light emission at 366 nm. Fermentation of the lactose causes an acidification of the medium and consequent precipitation of the biliary salts followed by absorption of the neutral red. The coliform organisms cultivate with pink-red colonies surrounded by a zone of precipitation while the organisms which do not ferment lactose cultivate with colorless colonies. The biliary salts and the crystal violet inhibit the growth of the Gram-positive bacteria.

## CULTURAL CHARACTERISTICS AFTER 18-24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Colonies colour	Fluore scence
<i>Ent. aerogenes</i>	good to excellent	pink-red without precipitation halo	-
<i>E. coli</i>	good to excellent	pink-red with precipitation halo	+
<i>Salm. enteritidis</i>	good to excellent	colourless	-
<i>Sh. dysenteriae</i>	good	colourless	-
<i>Staph. aureus</i>	inhibited	-	-
<i>Strep. faecalis</i>	inhibited	-	-

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