SELENITE BROTH

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pH = 7,0 ± 0,1

PREPARATION

Suspend 23 gr of powder in 1 litre of distilled or deionized water. Heat until completely dissolved. Dispense into sterile Petri dishes. Do not sterilize in the autoclave and avoid overheating.

USE

SELENITE BROTH is a enrichment medium for the isolation of *Salmonella* from faeces, urine or other pathologic material. In the first 6-12 hours of incubation the sodium selenite in the medium slows down the proliferation of the greater part of the intestinal commensals. To test faeces inoculate 1 gr of specimen in 8-10 ml of medium. For urines the medium must be used at double the concentration divided in 5-7 ml parts and inoculated with an equal volume of specimen.

CULTURAL CHARACTERISTICS AFTER 18-24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth
<i>Escherichia coli</i>	partially to completely inhibition
<i>Salmonella choleraesuis</i>	good to excellent
<i>Salmonella typhi</i>	good to excellent
<i>Salmonella typhimurium</i>	good to excellent

FORMULA (grams per litre)

Tryptone.....	5
Lactose.....	4
Sodium Phosphate.....	10
Sodium Selenite.....	4

pH = 7,0 ± 0,1

PRODOTTO PRODUCT	CODICE CODE	CONFEZIONE PACKAGING
SELENITE BROTH	610145	500 g. DCM
SELENITE BROTH	620145	100 g. DCM
SELENITE BROTH	6101455	5 kg. DCM
SELENITE BROTH	20110	10 PV
SELENITE BROTH	402050	6 FL x 100ml
SELENITE BROTH	412050	6 FL x 200ml

DCM = TERRENI DI COLTURA
DISIDRATATI / DEHYDRATED
CULTURE MEDIA

PS = PIASTRE PETRI / PETRI
DISHES

PV = PROVETTE / TUBES

FL = FLACONI / BOTTLES

BIBLIOGRAPHY

- Leifson E. (1963) - *New Selenite Enrichment Media for the isolation of Typhoid and Paratyphoid (Salmonella) Bacilli*. *Am. J. Hyg.* 24, 423-432.
- APHA (1963) - *Diagnostic Procedures and Reagents*, 4th Ed.