SIMMONS CITRATE AGAR

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PREPARATION

Suspend 24,2 gr of powder in 1 litre of distilled or deionized water. Heat until completely dissolved. Dispense in tubes. Sterilize in the autoclave at 121°C for 15 minutes. Allow the tubes to solidify in slant position.

USE

SIMMONS CITRATE AGAR is a medium used for differentiation and identification of Enterobacteria on the basis of their ability to utilize citrate as the sole source of carbon. Sodium citrate gets degraded to pyruvic acid by means of an enzyme system activated by bivalent cation supplied by magnesium phosphate. The medium also contains monoammonium phosphate as sole nitrogen source. The microorganisms which utilize citrate as sole source of carbon also utilize ammonium salts. Their degradation causes an alkalinity of the medium with color change of the indicator to blue. The citrate positive bacteria cultivate with color change of the indicator to blue while the citrate negative bacteria do not grow.

CULTURAL CHARACTERISTICS AFTER 18-24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Medium colour
<i>Enterobacter aerogenes</i>	good	blue
<i>Escherichia coli</i>	inhibited	green
<i>Salmonella enteritidis</i>	good	blue
<i>Salmonella typhimurium</i>	good	blue
<i>Salmonella typhi</i>	good	green
<i>Shigella dysenteriae</i>	inhibited	green

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BIBLIOGRAPHY

1. Edwards, P.R. and Ewing, W.H. (1972). *Identification of enterobacteriaceae*. 3rd Edition, Minneapolis: Burgess Publishing Company.
2. Mac Faddin J.F. (1976). *Biochemical tests for identification of medical bacteria*. Baltimore: The Williams and Wilkins Company.