T.S.I. AGAR

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PREPARATION

Suspend 64,5 gr of powder in 1 litre of distilled or deionized water. Heat until completely dissolved. Dispense in tubes and sterilize in the autoclave at 121°C for 15 minutes. Cool in inclined position so as to obtain a base at least 2 cm deep.

USE

T.S.I. AGAR is a medium used for the differentiation of enterobacteria according to lactose, glucose and sucrose fermentation and to hydrogen sulphide and carbon dioxide production. Fermentation of the sugar, resulting in acid production is detected by the phenol red indicator, the color changes to yellow and red for alkalinization. Fermentation of lactose and sucrose determines an acid reaction on the surface and in depth. Lack of sugars fermentation determines an alkaline reaction both in depth and on the surface. Sodium thiosulphate is reduced to hydrogen sulphide which then reacts with an iron salt yielding the typical black iron sulphide. Gas production is determined by the formation of bubbles up to a more or less severe fragmentation of the agar.

REACTION'S TABLE

Microorganism	Slant	Bottom	H ₂ S
<i>Escherichia</i>	A (K)	A/G or A	-
<i>Shigella</i>	K	A	-
<i>S.typhi</i>	K	A	±
<i>Altre Salmonella</i>	K	A/G	+
<i>Arizona</i>	A (K)	A/G	+
<i>Arizona</i>	A (K)	A/G	+
<i>Citrobacter</i>	K	A/G	+
<i>Edwardsiella</i>	A	A/G	-
<i>Klebsiella</i>	A	A/G	-
<i>Enterobacter</i>	K	A/G	-
<i>E. Hafniae</i>	K or A	A	-
<i>Serratia</i>	A (K)	A/G	+
<i>Proteus vulgaris</i>	K (A)	A/G	+
<i>P. Mirabilis</i>	K	A/G or A	-
<i>P. morganii</i>	K	A	-
<i>P. rettgeri</i>	K	A/G or A	-

K = alkaline reaction

A = acid reaction

A/G = acid reaction with gas production

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CULTURAL CHARACTERISTICS AFTER 24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Slant	Bottom	Gas	H ₂ S
<i>E. coli</i>	good	yellow	yellow	+	-
<i>C. freundii</i>	good	yellow	yellow	+	+
<i>P. vulgaris</i>	good	yellow	yellow	+	+
<i>S. enteritidis</i>	good	red	yellow	+	+
<i>S. flexneri</i>	good	red	yellow	-	-
<i>P. aeruginosa</i>	good	red	red	-	-

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BIBLIOGRAPHY

1. Edwards P.R. & Ewing W.H. (1972). *Identification of enterobacteriaceae*, 3rd edition, Minneapolis: Burgess Publishing Company.
2. Mac Faddin J.F. (1976). *Biochemical tests for identification of medical bacteria*. Baltimore: The Williams & Wilkins Company.