VIOLET RED BILE MUG AGAR

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PREPARATION

Suspend 41,6 gr of powder in 1 litre of distilled or deionized water. Heat to boiling while shaking. Cool to 50°C and transfer into sterile Petri dishes. DO NOT STERILIZE IN THE AUTOCLAVE. AVOID OVERHEATING.

USE

VIOLET RED BILE MUG AGAR is a selective and differential medium for the isolation and enumeration of coliform bacteria in foods, milk and dairy products. The MUG (4-methylumbelliferil-β-D-glucuronide) allows detection of *Escherichia coli* directly in the Petri dish. In fact MUG, which is colorless is hydrolyzed by the enzyme β-glucuronidase in the *Escherichia coli* while releasing the compound 4-methylumbelliferone, highly fluorescent if observed under a UV 1 amp with light emission at 366 nm. Biliary salts and crystal violet inhibits the growth of the Gram-positive bacteria. The neutral red color allows detection of lactose-fermenting bacteria from the non-fermenting ones. Lactose fermentation causes acidification of the medium with consequent color change of the indicator towards red-purple and with salts precipitation.

CULTURAL CHARACTERISTICS AFTER 24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Colonies colour	Fluor.
<i>Ent. aerogenes</i>	good to excellent	pink w/o precipitation halo	-
<i>E. coli</i>	good to excellent	pink-violet w precipitation halo	+
<i>Staph. aureus</i>	inhibited	-	
<i>Strept. faecalis</i>	inhibited	-	

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BIBLIOGRAPHY

1. APHA (1985). *Standars methods for the examination of water and wastewater*, 16th edition, Washington, D.C.
2. Robinson, B.J. (1984). *App. Environ. Microbiol* 48, 285-288.