VOGEL JOHNSON AGAR

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PREPARATION

Suspend 60 gr of powder in 1 litre of distilled or deionized water. Heat until completely dissolved. Sterilize in the autoclave at 121°C for 10 minutes. Cool to 50°C and aseptically add 20 ml of a 1% potassium tellurite sterile solution. Mix well and dispense into sterile Petri dishes.

USE

VOGEL JOHNSON AGAR is a selective medium used for the isolation of *Staphylococcus aureus* from food and clinical specimens. Mannitol and red phenol allow to differentiate the mannitol fermented staphylococci. Mannitol fermentation is indicated by a yellow zone around the colonies. The medium's selectivity is due to glycine, lithium chloride and potassium tellurite which inhibit the growth of Gram-positive and Gram-negative bacteria except for *Staphylococcus aureus*. When potassium tellurite is reduced to tellurium it forms a black precipitation inside the colonies. Members of the *Proteus* strains cultivate after 18 hours with black colonies and with a modification of the medium towards brown. After 48 hours of incubation an alkaline inversion of the pH of the medium occurs resulting in a purple color.

CULTURAL CHARACTERISTICS AFTER 24-48 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Colonies colour
<i>Escherichia coli</i>	inhibited	
<i>Proteus mirabilis</i>	poor	black
<i>Staphylococcus aureus</i>	poor	black with yellow halo
<i>Staph. epidermidis</i>	good	translucent or black

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BIBLIOGRAPHY

1. APHA (1976). *Compendium of methods for the microbiological examination of foods.*
2. Vogel Johnson M (1960). *A modification of the tellurite glycine medium for use in the identification of Staphylococcus aureus. Public Health Lab., 18, 131.*