ENDO AGAR LES

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PREPARATION

Suspend 51 gr. of powder in 1 litre of distilled or deionized water. Add 20 ml of ethyl alcohol 95°. Bring to boil while gently shaking. Cool to 50°C. Dispense into sterile petri dishes. Do not sterilize in the autoclave and avoid overheating

USE

ENDO AGAR LES is a medium recommended for the enumeration of coliform organisms in water by the membrane filter technique using LAURYL SULPHATE BROTH as a preliminary enrichment. Sodium desoxycholate and sodium lauryl sulphate inhibit the growth of Gram-positive bacteria. The coliform organisms cultivate with red colonies having a distinctive metallic appearance.

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CULTURAL CHARACTERISTICS AFTER 24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Colonies colour
<i>Enterobacter aerogenes</i>	good to excellent	red
<i>Escherichia coli</i>	good to excellent	metallic red
<i>Salmonella typhi</i>	good to excellent	colourless
<i>Shigella sonnei</i>	good to excellent	colourless

FORMULA (grammi per litro)

Estratto di Lievito.....	1.2
Tryptone.....	3.7
Peptone.....	2.5
Tryptospecial.....	7.5
Lattosio.....	9.4
Potassio Fosfato dibasico.....	3.3
Potassio Fosfato monobasico.....	1
Sodio Cloruro.....	3.7
Sodio Desossicolato.....	0.1
Sodio Laurilsolfato.....	0.05
Sodio Solfito.....	1.6
Fucsina Basica.....	0.8
Agar.....	15
pH = 7.2 ± 0.2	

FORMULA (grams per litre)

Yeast Extract.....	1.2
Tryptone.....	3.7
Peptone.....	2.5
Tryptospecial.....	7.5
Lactose.....	9.4
Dipotassium Phosphate.....	3.3
Monopotassium Phosphate.....	1
Sodium Chloride.....	3.7
Sodium Desoxycholate.....	0.1
Sodium Lauryl Sulphate.....	0.05
Sodium Sulphite.....	1.6
Basic Fuchsin.....	0.8
Agar.....	15
pH = 7.2 ± 0.2	

PRODOTTO PRODUCT	CODICE CODE	CONFEZIONE PACKAGING
ENDO AGAR LES	611011	500 g. DCM
ENDO AGAR LES	621011	100 g. DCM
ENDO AGAR LES	11053	20 PS
ENDO AGAR LES	11053*	100 PS
ENDO AGAR LES	16343	10 PS

DCM = TERRENI DI COLTURA
DISIDRATATI / DEHYDRATED
CULTURE MEDIA

PS = PIASTRE PETRI / PETRI
DISHES

PV = PROVETTE / TUBES

FL = FLACONI / BOTTLES

BIBLIOGRAPHY

1. McCarthy, De Laneg and Grasso (1961). *Water and Sewage Works*, 108, 238.
2. APHA (1985). *Standards Methods for the Examination of Water and Wastewater*, 16th Ed.