# HEKTOEN ENTERIC AGAR

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pH = 7.4 ± 0.2

## PREPARATION

Suspend 78 gr. of powder in 1 litre of distilled or deionized water. Bring to boil while gently shaking. Cool to 50°C. Dispense into sterile Petri dishes. Do not autoclave. Avoid overheating.

## USE

HEKTOEN ENTERIC AGAR is a differential and selective medium for the isolation of Gram-negative intestinal flora. Biliary salts inhibit the growth of the Gram-positive bacteria. The differentiation of the intestinal flora is based on the following biochemical tests:

- 1) fermentation of lactose, sucrose and salicin decreases pH of the medium and precipitates biliary salts. The fermenting colonies take on a salmon color while the non fermenting ones cultivate with green-blue colonies;
- 2) the production of hydrogen sulfate initiated by the sodium thiosulfate induces the precipitation of ferric sulphide causing colonies with a black centre.

## CULTURAL CHARACTERISTICS AFTER 24 HOURS OF INCUBATION AT 37°C

Microorganism	Growth	Colonies colour
<i>Shigella flexneri</i>	excellent	green-blue
<i>Salmonella enteritidis</i>	excellent	green-blue with or without H <sub>2</sub> S production
<i>Salmonella arizona</i>	excellent	green-blue with or without H <sub>2</sub> S production
<i>Escherichia coli</i>	good	yellow
<i>Proteus spp</i>	excellent	green-blue with or without H <sub>2</sub> S production

## FORMULA (grams per litre)

Peptone.....	12
Yeast Extract.....	4
Bile Salts n.3.....	9
Lactose.....	12
Saccharose.....	12
Salicin.....	2
Sodium Chloride.....	5
Sodium Thiosulphate.....	5
Ferric Ammonium Citrate.....	1.5
Brom Thymol Blue.....	0.065
Acid Fuchsin.....	0.1
Agar.....	15

pH = 7.4 ± 0.2

<b>PRODOTTO PRODUCT</b>	<b>CODICE CODE</b>	<b>CONFEZIONE PACKAGING</b>
HEKTOEN ENTERIC AGAR	610021	500 g. DCM
HEKTOEN ENTERIC AGAR	620021	100 g. DCM
HEKTOEN ENTERIC AGAR	6100215	5 kg. DCM
HEKTOEN ENTERIC AGAR	10043	20 PS
HEKTOEN ENTERIC AGAR	10043*	100 PS
HEKTOEN ENTERIC AGAR	402230	6 FL x 100ml
HEKTOEN ENTERIC AGAR	412230	6 FL x 200ml

DCM = TERRENI DI COLTURA PS = PIASTRE PETRI / *PETRI* PV = PROVETTE / *TUBES* FL = FLACONI / *BOTTLES*  
 DISIDRATATI / *DEHYDRATED* DISHES  
 COLTURE MEDIA

**BIBLIOGRAPHY**

1. King, S.& Metzger, W.I. (1968). A New Plating Medium for the Isolation of Enteric Pathogens. *Hektoen Enteric Jagar, Appl. Microbiol.*, 16, 577-578.
2. King, S.& Metzger, W.I. (1968). A New Plating Medium for the Isolation of Enteric Pathogens. Comparison of *Hektoen Enteric Agar With SS and EMB Agar, Appl. Microbiol*, 16, 579-581