

Simboli usati per i Diagnostici in Vitro

KEY GUIDE TO SYMBOLS SIGNIFICATION DES SYMBOLES ERLÄUTERUNG DER SYMBOLE CLAVE DE LOS SÍMBOLOS UTILIZADOS GUIDA AI SIMBOLI.




Use by
A utiliser avant le
Verw. bis:
Utilizar antes de
Usar entro.



Lot Lot
Ch.-B.:
Lote
Lotto.


REF

Catalogue number
No. dans le catalogue
Katalognummer
No de catálogo
No. di catalogo.



For *in vitro* diagnostic use
Pour usage diagnostique *in-vitro*
In-vitro diagnostikum
Para uso diagnóstico *in vitro*
Per uso diagnostico *in vitro*.

Consult , accompanying documents
Voir documents cijoins
Siehe beigefügte, Dokumente
Véanse los documentos djuntos
Vedi documenti allegati



8°C

Store at 2-8°C
Conserver à 2-8°C
Lagern bei 2-8°C
Conservar a 2-8°C
Conservare a 2-8°C



Biological risks
Risques biologique
Biogefährdung
Riesgo biológico
Rischio biologico



Buffer Solution
Solution Tampon
Pufferlösung
Solución Tampón
Soluzione Tampone

Recon.

Reconstitute with
Reconstituer avec
Lösen in
Reconstituya con
Ricostituire con



Positive Control
Contrôl positif
Positive kontrolle
Control positivo
Controllo positivo



Negative Control
Contrôle négatif
Negative kontrolle
Control negativo
Controllo negativo