

	ENTEROSYSTEM SYSTEM – 18 R - REAGENT Reagenti biochimici per Enterosystem 18 R <i>Biochemical reagents for Enterosystems 18 R</i>	CODE 80252
--	---	-------------------

REAGENTE / REAGENT	CONTENUTO / CONTENT
ALPHA NAPHTHOL Test Voges-Proskauer <i>Voges Proskauer test</i>	<i>Alpha naphthol</i>0.6 g <i>Ethyl alcohol</i> 10 ml
KOVAC'S Test indolo <i>Indole test.</i>	<i>p-Dimethylaminobenzaldeide</i>0.5 g <i>Amilic alcohol</i>7.5 ml <i>Hydrochloric acid 37%</i>2.5 ml
VASELIN OIL Reagente per reazioni anaerobiche <i>Reagent for anaerobic reactions</i>	<i>Sterile vaselin oil</i>40 ml

CONTENUTO (100/200 tests) 3 Flacons x 10 ml 2 Flacons x 40 ml CONSERVAZIONE: 5-12°C VALIDITA': 1 anno	CONTENT (100/200 tests) 3 Bottles x 10 ml 2 Bottles x 40 ml STORE AT: 5-12°C SHELF LIFE: 1 year
---	--

--	--