

PKL PPC 175

Fully-auto Blood
Coagulation Analyzer

- King of emergency
- Unique in the world
- Accurate, convenient



PKL PPC 175

Fully-auto Blood Coagulation Analyzer

Specifications

Test method: Coagulation method, immunoturbidimetry optional

Test item: APTT, PT, TT, FIB and blood coagulation factors

Coagulation-related special items D-Dimer and FDP optional

PT-derived FIB function

Test speed: 100 T/H

Sample position: 6 sample positions

Tilt Reagent position: 11 reagent positions

Test position: 4 positions

Cleaning position: 1 position

Sample probe: 1 probe

Reagent system: built-in card reader

Functions:

Reagent position interchange function

Reagent position refrigeration function

Sample probe constant heating function

Sample liquid level detection function

Sample and reagent shortage alarm function

Four directions ventilation cooling function

LED work lighting function

Features

Automatic and manual dual-model test capability

Tilt reagent position to save reagent

Display coagulation curve of dynamic test

Automatic and manual calibration function

Batch test and emergency priority insertion test function

Abnormal test results alarm and automatic retest function

Background scan of optical coagulation method, remove jaundice, high fat background interference function

Cleaning solution shortage alarm and waste liquid overflow alarm function

PARAMEDICAL srl

Via Irno snc,

84098 Pontecagnano Faiano (SA) - ITALY

Vat.: IT03117920656

Tel.+39 089.385027

Fax+39 089.3854479

amministrazione@paramedical.it

www.paramedical.it

