

Chemistry Analyzer

PKL PPC 115 Chemistry Analyzer

- High performance
- Easy touch screen operation
- Convenient maintenance design
- Support barcode reader
- With 8 inner incubation positions

Specifications

- Analysis method: Kinetic, end point, two point, dual wavelength, ABS
- Light source: halogen lamp, 6V/10W Wavelength (in nm): 340, 405, 492, 510, 546, 578, 630, 700
- Absorbance range: -0.3~4.0 Abs
- Wavelength accuracy: ± 2 nm
- Spectrum band: ± 10 nm




- Flowcell: 10mm quartz cell, 32 μ L volume
- Temperature control: room temperature, 37C (± 0.1 C variable)
- Sampling volume: 100 μ L ~ 9999 μ L adjustable, recommended: 500 μ L
- Cross contamination: $\leq 1.0\%$
- Stability: ≤ 0.005 Abs/h
- Memory: Up to 300 items and 30000 test Results

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

- Display: Large LCD touch screen display
- Printer: Built-in thermal printer
- Interface: Parallel port (For external printer) ,RJ-45,3USB,SDID,PS-2
- Ambient: Room temperature: 10 C~32C
- Humidity: ≤40%~85%
- Power: AC (110-220)V, 50Hz
- Dimension:362.5 x 341x 165mm (length * width * height)



Features

- High-quality grating & close optical system
- Easy touch screen operation with 7 inch colour screen
- Convenient maintenance design
- Large memory capacity
- With 8 inner incubation positions
- Support barcode reader

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

Via Antonio Amato , 24
84131 SALERNO - ITALY
Vat.: IT03117920656
Ph + 39 089 385027 amministrazione@paramedical.it www.paramedical.it

