

## Electrolyte Analyzer

### **PKL PPC 182** Electrolyte Analyzer

#### Features

- K, Na, Cl, Ca, Li, CO<sub>2</sub>, pH, Mg, Ag
- 7" Color touch screen, easy to operate
- High accurate and long life electrode and TCO<sub>2</sub> sensor
- Programmable multi-format print-out
- Reagent pack, real time monitoring of reagent residual volume
- Up to 50,000 test results can be stored
- RS-232 port, supporting bar code reader and USB interface
- Sleep mode to reduce reagent consumption

#### Options:

- Auto Sampler
- Self-made ISE, maintenance free and long life span.
- Original reagent, proven quality and one source convenience



**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**





**PKL**<sup>®</sup>  
POKLER ITALIA

**Parameter configuration (models):**

- A** K, Na, Cl
- B** K, Na, Cl, TCO2
- C** K, Na, Cl, iCa, nCa, TCa, pH
- D** K, Na, Cl, iCa, nCa, TCa, pH, TCO2, AG
- F** K, Na, Cl, Li
- G** K, Na, Cl, Li, TCO2
- H** K, Na, Cl, iCa, nCa, TCa, pH, Li
- I** K, Na, Cl, iCa, nCa, TCa, pH, Li, TCO2, AG
- J** K, Na, Cl, iCa, nCa, TCa, pH, Mg, TCO2, AG
- K** K, Na, Cl, iCa, nCa, TCa, pH, Li, Mg, TCO2, AG

**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**

Via Antonio Amato , 24  
84131 SALERNO - ITALY  
Vat.: IT03117920656  
Ph + 39 089 385027 [amministrazione@paramedical.it](mailto:amministrazione@paramedical.it) [www.paramedical.it](http://www.paramedical.it)

