

Hematocrit Centrifuge

PKL PPC 512 Microfuge-24 Centrifuge

Features

- This unit is designed for flexibility and ease of use in the laboratory, hospital and blood bank.
- Brushless frequency motor for model PPC 512, free maintenance, no powder pollution, quick in speed up and down.
- Emergency lid opening switch, rubber fastenings, small vibration and vibration-free operation.
- There are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.

Micro computer control system, LCD display the RCF time and speed.



PPC 512 Technical Data:

- Max speed: 12000rpm
- Speed accuracy: ± 10 rpm
- Max RCF: 15800 \times g
- Noise: ≤ 55 dB(A)
- Max volume: 12 capillaries and 24 capillaries rotors for your choice.
- Power supply: AC220V 50HZ
- Timer: 0~99min
- Safety shut-off switch in lid, compact design, super speed and imbalance protection.
- The centrifuge body is made of high-quality steel, very safe and reliable.
- Dimension(LxWxH): 355 \times 270 \times 205mm
- Weight: 13KG
- Package: Carton box

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

Via Antonio Amato , 24
84131 SALERNO - ITALY
Vat.: IT03117920656
Ph + 39 089 385027 amministrazione@paramedical.it www.paramedical.it

