

Class I Biosafety Cabinets

PKL PPC BIO 100/101

Class I: The Class I biological safety cabinet is an open-front negative pressure cabinet. The Class I biological safety cabinet will provide personnel and environmental protection, but not product protection. All of the air from the cabinet is sent through a HEPA filter, either into the laboratory or to the outside. Optional: Base stand



<u>MODEL</u>	<u>Internal Dims</u>	<u>External Dims</u>	<u>Package Size</u>
PPCBIO100	540x385x440	550x395x730	755x600x950
PPCBIO101	680x540x500	700x550x900	1000x700x1200

PKL PPC BIO 102

Class I: The Class I biological safety cabinet is an open-front negative pressure cabinet. The Class I biological safety cabinet will provide personnel and environmental protection, but not product protection. All of the air from the cabinet is sent through a HEPA filter, either into the laboratory or to the outside. Optional: Base stand



<u>MODEL</u>	<u>Internal Dims (mm)</u>	<u>External Dims (mm)</u>	<u>Package Size (mm)</u>
PPCBIO102	540x385x440	550x395x730	755x600x950

PKL PPC BIO 103

Class I: The Class I biological safety cabinet is an open-front negative pressure cabinet. The Class I biological safety cabinet will provide personnel and environmental protection, but not product protection. All of the air from the cabinet is sent through a HEPA filter, either into the laboratory or to the outside. Optional: Base stand



<u>MODEL</u>	<u>Internal Dims (mm)</u>	<u>External Dims(mm)</u>	<u>Package Size(mm)</u>
PPC BIO 103	968x695x630	1100x695x1924	1270x880x2230

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016



PKL[®]
POKLER ITALIA

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

Via Antonio Amato , 24
84131 SALERNO - ITALY
Vat.: IT03117920656
Ph + 39 089 385027 amministrazione@paramedical.it www.paramedical.it

