



VETERINARY PRODUCTS

PKL PPC 1200H VET

5-Diff Hematology Analyzer

- 28 parameters including RET%, RET#, IRF
- Laser light multi-dimensional cell classification
- Double mode for WBC counting (Impedance Count and Optical Count)
- Resistant RBC mode
- Up to 60T/H
- Internal barcode reader



PKL PCC 1040H VET

5-Diff Hematology Analyzer

- 28 parameters including RET%, RET#, IRF
- 60T7H
- Laser light multi-dimensional cell classification
- Resistant RBC mode
- Support both whole blood and capillary blood samples
- Only three reagents for 25 parameters



PKL PCC 1050H VET

5-Diff Hematology Analyzer

- 28 parameters including RET%, RET#, IRF
- 60T7H
- Laser light multi-dimensional cell classification
- Resistant RBC mode
- Support both whole blood and capillary blood samples
- Only three reagents for 25 parameters



PKL PCC 610H VET

3-Diff Hematology Analyzers

- 16 pre-determined mode settings
- 60T/H, 21 parameters and 3 histograms
- 10.4 inch color LCD display & Linux operation system
- Support USB data backup & system upgrade
- Support LIS and HIS with HL7 protocol
- Touch screen optional



ELPHO VET

Automatic Electrophoresis System on Acetate Cellulose

- Analysis type: serumprotein, hemoglobin
- 1 to 16 samples micro
- 16 samples version: in 60 min; the first 8 samples in 40 min
- 8 independent channels
- Chambers can be both extractable in order to be filled.
- Emptied and washed; or fixed with a total fluidic automation
- Computer inside, USB, LAN.
Remote connection.



ELPHO AGAROSE VET

Automatic Electrophoresis System on Agarose Gel

- Serumproteins, Hemoglobins, Urine Proteins
- 1 to 26 samples: 2 lines of 13 samples (total 26 samples holes); 26 samples in 45 min
- 13 ultrabright led
- Load, recover and discharge all reagent with fluidic full Automatic
- Temperature Migration controller with Peltier and transistor for dessication of gel
- In the internal small basins with time infrared probes;
- in the external tanks with immersion probes.
- Windows 10 or Windows 7
- PC & Keyboard ; mouse. (external PC)
- Dedicated Software with methodic changeable management



PKL PPC 100 VET

Urine Analyzer

- Up to 120T/H and
14 parameters available
- FDA approved
- Portable
- Barcode reader optional
- Urin-10K, Urin-11K and
Urin-14K control available



PKL PPC 110VET

Chemistry Analyzers

- Easy visual reading
- Much higher reliability with
Ascorbic Acid parameter



PKL PCC 115VET

Chemistry Analyzer

- Friendly operation with 7 inch colourful touch screen
- Sophisticated Software, large memory capacity
- High-quality filters & close optical system
- Convenient maintenance design



PKL PPC 125VET

Automatic Chemistry Analyzer

- Carry-on baggage size
- End point $cv \leq 0.5\%$; Kinetic $cv \leq 1\%$, two points $cv \leq 1.5\%$
- Lower failure rate; easy maintenance
- 9 Wavelengths
- Barcode reader (optional)
- Water cycle cooling
- Auto warm water washing
- Immunology items
- Powerful software
- Support lis interface
- 200 tests/hour



PKL PPC 200 VET

Automatic Chemistry Analyzer

- Random Access, constant 200 tests/hour
- 24 hour non-stop cooling system to ensure reagent at 2-8°C
- Durable ceramic syringes to ensure accuracy & precision
- High accurate optical system
- Collision protection in both vertical and horizontal directions, stop & alarm automatically once touching barrier, not affect former tests
- Automatic eligible cuvettes detection & selection
- Select best test point by reaction curve, create new factor automatically
- Support LIS interface



PKL PCC 9000 VET

Integrated Diagnostic System

- General clinic information technology solution
- Wide angle detection
- Detachable disinfection
- Combination design
- Support bluetooth

