

PKL PPC 175/175 Plus

Fully-auto Blood
Coagulation Analyzer



- King of emergency
- Unique in the world
- Accurate, convenient



PKL PPC 175/175 Plus

Fully-auto Blood Coagulation Analyzer

Specifications

Test method: Photoelectric colorimetry

Test item: APTT, PT, TT, FIB

D-Dimer optional (PKL PPC 175/175 Plus)

PT-derived FIB function

Test speed: 100 T/H

Sample position: 5 sample positions

Tilt Reagent position: 6 reagent positions

Test position: 4 positions

Cleaning position: 1 position

Sample probe: 1 probe

Reagent system: external card reader

Functions:

Reagent position refrigeration function optional (PKL PPC 175/175 Plus)

Sample probe constant heating function

Sample liquid level detection function

Sample and reagent shortage alarm function

Four directions ventilation cooling function

LED work lighting function

Features

Automatic and manual dual-model test capability

Tilt reagent position to save reagent

Display coagulation curve of dynamic test

Automatic and manual calibration function

Batch test and emergency priority insertion test function

Abnormal test results alarm and automatic retest function

Background scan of optical coagulation method, remove jaundice, high fat background interference function

Cleaning solution shortage alarm and waste liquid overflow alarm function

PARAMEDICAL srl

Via Antonio Amato, 24 - 84131 Salerno (Italy)

Vat.: IT03117920656

Tel.+39 089.385027

Fax+39 089.3854479

amministrazione@paramedical.it

www.paramedical.it

