

Micro High Speed Centrifuge

PKL PPC 503 Biofuge PCR

Features

- 1.Brushless motor, no pollution, free-maintenance.
- 2.Microprocessor control, LED display which indicator the speed, time, RCF inoperation,10 kind sof brake setting, operate simply.
- 3.Electric lid lock, superspeed, imbalance protection. The centrifuge body is made of high quality steel, stainless steel chamber, safe and reliable.
- 4.Rotor is connected to spindle by specialized taper sleeve ,loading simple and quick, no direction.
- 5.3 tiers protection steel cover and get the ideal centrifugation result



PPC 503 Technical Data:

- | | |
|-----------------------|--------------------------|
| • Max Speed | 16000 r/min |
| • Max Capacity | 10 x 5 mL |
| • Max RCF | 19040 x g |
| • Timer | 1 - 99h 59 min |
| • Acc / Dec | 10 kinds |
| • Noise | ≤ 55 dB (A) |
| • Power | 110/220 V / 60 Hz / 2A |
| • Dimensions | 355 mm x 270 mm x 205 mm |
| • Net Weight | 14 Kg |

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

Rotor detail information:

Rotor No.	Max Speed(rpm)	Max Volume(ml)	Max.RCF (×g)
Angle Rotor	16000	40×0,2ml	19040
Angle Rotor	16000	24×0,5ml	18480
Angle Rotor	16000	12×1,5/2ml	16260
Angle Rotor	16000	10×5ml	17880
Angle Rotor	14000	20×1,5/2ml	15580
PCR Rotor	13000	4x8 PCR	10410



NO31003
Angle Rotor
12×1.5/2ml
16000rpm



NO61006
PCR Rotor
4×8PCR
13000rpm

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016