

5 Digital 3 D Shaker

PKL PPC 3D - 5

5 Digital 3 D Shaker

Features

3D shakers are very much essential instruments for mixing homogeneous mixing in flasks, culture flasks, Petri dishes and tubes. Ideal for e.g. the cultivation of cell cultures, distributing of cells and staining and de-staining of cells. Neuation has come up with 3kg and 5kg capacity 3D shakers with three-dimensional tumbling motion and adjustable speed for smooth mixing tasks. These Digital 3D shaker comes with PUSLE mode feature and FAV mode.record data logs



Specifications

Suitable for the mixing of upto 5kg Capacity
Tumbling rocking motion -7 deg. tilt
Speed range of 5 to 120 RPM for gentle to fast motion and timer function
Long life BDLC motor with Microprocessor control & digital display
Unique Run Pause Run mode, user can set 30 to 99 sec
Platform size - 307 x 267 mm - with spiked mat
User can set and save upto 9 user defined programs
IP 21 class protection with removable platform for ease of cleaning
Standard Accessory included: <u>Spiked Rubber Mat</u>
Spiked rubber mat (Standard)
(Optional) - Sticky Mat
Stacking/ double decker platform with 4 posts and spiked rubber mat

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016