

Overhead Stirrer Elite Series

**PKL PPC PPC Elite**  
Overhead Stirrer

**Overhead Stirrer Elite Series**

PKL PPC Elite Series

Features

- Brushless DC motor for quiet operation and maintenance free long operation
- Microprocessor controlled for maintenance of accuracy even with viscous liquid
- Viscosity up to 70,000 mPas & 100 L volume
- Adjustable speed ranging from 50 Up to 2000 RPM
- Long timer setting from 1min to 99 hours 59 mins & continuous mode
- With reversible mixing motion - CW/CCW
- Large digital display for all the parameters at one glance
- Uninterrupted operation through Smart control
- Safety cut for overload, over current & over temperature - safety cut



**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**



<b>Specifications</b>	<b>Stir OH 50 Elite</b>	<b>Stir OH 75 Elite</b>	<b>Stir OH 100 Elite</b>
Motor Type	Brushless DC	Brushless DC	Brushless DC
Capacity	50 L	75 L	100 L
Speed Range	50-2000 RPM	50-2000 RPM	50-1300 RPM
Torque	50 Ncm	75 Ncm	100 Ncm
Viscosity	30000 mPas	50000 mPas	70000 mPas 1 min to 99 hr 59 mins
Timer	1 min to 99 hr 59 mins & infinity mode	1 min to 99 hr 59 mins & infinity mode	1 min to 99 hr 59 mins & infinity mode
Direction of Rotation	Reversible - CW/CCW	Reversible - CW/CCW	Reversible - CW/CCW
Protection Class	IP 42	IP 42	IP 42
Dimension(Head)	102x 249 x 242mm	102x 249 x 242mm	102x 249 x 242mm
Weight	5.1 kgs	5.1 kgs	5.1 kgs



**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**

