

Smart Personal Micro Centrifuge

PKL PPC M08 Smart Centrifuge

Features

- ☐ 8 place closed rotor ensures better airflow, greater efficiency with quiet operation and low heat generation ensures stability for temperature sample
- ☐ Speed Accuracy +/- 100 RPM
- ☐ Small footprint saves valuable bench space - suitable for use in hoods and cold rooms
- ☐ Electronic brake for immediate stop on lid opening for user safety
- ☐ Saves time through quick runs by simply closing and opening the lid when switched On
- ☐ Carbon Dust insulation and safe drainage provision in design



Specifications

Rotor Capacity

Speed RPM

Speed RCF (g)

Digital Display

Electronic Safety Brake

Noise Level

Power Supply

Dimensions (WxDxH) in mm

Rubber suction pads

Rotors and Adaptors

(included in standard pack)

PKL PPC M08

8 Place Closed Rotor or 1.5/2.0 ml microtubes

6000

2000

Yes

Yes

< 55 db

Universal (95 to 230V, 50/60 Hz)

162 x 157 x 115 mm

Yes

8 Slot x 1.5/2.0 ml Closed Rotor

PCR Strip Rotor (2 x 8 x 0.2 ml) Reduction

Adaptors for 0.2 / 0.4 ml Microtube

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

Via Antonio Amato , 24
84131 SALERNO - ITALY

Vat.: IT03117920656

Ph + 39 089 385027 amministrazione@paramedical.it www.paramedical.it

