

## Overhead Stirrer

# PKL PPC PPC OHJr Overhead Stirrer

### Features

- Brushless DC motor for quiet operation and maintenance free long operation
- Microprocessor controlled for maintenance of accuracy even with viscous liquid
- Stirs viscosity up to 1,000 mPas or 10 L volume capacity (H<sub>2</sub>O)
- Adjustable speed ranging from 50 Up to 2200 RPM
- Safety cut for overload
- Plug & play design for easy usage



### Specifications

Motor Type  
Capacity(H<sub>2</sub>O)  
Speed Range  
Torque  
Viscosity  
Protection Class  
Weight

### Stirrer OH Jr.

Brushless DC  
10 L  
50-2200 RPM  
10 Ncm  
1000 mPas  
IP 42  
1.5 kgs

**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**