

## Magnetic Stirrer Penta

# PKL PPC Penta

## Magnetic Stirrer

### Features

- Brushless DC motor for quiet operation and maintenance free long life
- Corrosion resistant SS top surface for better chemical resistance
- Microprocessor controlled stirrer with variable speed and time setting and along with last run memory feature
- Adjustable speed from 100 to 2200 RPM in steps of 10 RPM
- Long timer ranges up to 99 hours and 59 minutes & infinite mode
- Programmable Pulse mode for homogenous stirring in clockwise and counter-clockwise direction
- Low aspect ratio (just 45 mm thick)
- Software driven gradual acceleration
- Universal power supply - small wall mount adaptor
- IP-42 compliant design



### Specifications

Motor type  
Capacity (water)  
Speed setting  
Timer  
Digital display  
Pulse mode  
Last run memory  
Dimensions (W x D x H)  
Protection class  
Weight

### *Magnetic Stirrer Penta*

Brushless DC  
10 litre (2 litre per station)  
100 to 2200 RPM (at every 10 RPM)  
Up to 99 hours & 59 minutes & infinite mode  
Yes  
Yes programmable (30 to 99 sec)  
Yes  
594x332x45 mm  
IP 42  
8 kg

**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**