

# Plate Shaker 4 Plates

## **PKL PPC SH4** 4 Plate Shaker

### Features

- 3mm orbital motion for thorough mixing in microplates
- Adjustable speed 200 to 1200 RPM
- Rubber O-ring provides firm grip for plate holding
- Easy access to all 4 plates for loading and unloading
- Programmable Pulse mode for vigorous shaking in clockwise and anticlockwise direction
- Last run memory feature
- Compatible to all standard skirted microplates
- Brushless DC motor for quiet operation and maintenance free long life
- Easy to clean and maintain



### Specifications

Motor Type	<b>SH4</b> BLDC motor
Orbital Diameter	3 mm
Number of Plates	4
Speed Setting	200 to 1200 RPM
Timer setting	1 to 999 mins & infinite mode
Pulse Mode	Yes, programmable (30 to 90 secs)
Maximum volume	4x384x60 µl
Weight	5.4 kg
Inbuilt Counterbalance	Yes
Dimension (W X D X H)	343 x 195 x 88 mm
Digital display	24V
Included Attachment	Acc. All standard Skirted Microplates

**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**

