

Vortex Mixer

PKL PPC VMA1 Vortex Mixer

Features

- Variable speed selection - delivers upto 2800 RPM
- 3 way switch (Touch, Continuous & Stand-by mode)
- Gentle press for smooth Operation
- Brushless AC motor fitted with ball bearing
- Orbital diameter of 4 mm
- Visual indication marker for max speed with attachment
- IP 21 class protection - suitable for use in hoods
- Optional attachments for low speed shaking with Microplate, Microtubes & Flask






Specifications

	VMA1
Max Speed	2800 RPM
Orbital diameter	4 mm
Motor type	Brushless AC
Max Load Capacity	500 gm
Speed accuracy	± 5%
Different modes Touch, Stand-by & Continuous mode	Touch, Stand-by & Continuous mode
Dimensions	205 x 136 x 138.5 mm
Power Supply	220 - 240VAC, 50Hz
Ambient temperature	5 to 40°C
Weight	3 kg
Optional Attachments (not included)	Universal platform for attaching either microplate or microtube or flask Attachments

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016



Optional Accessories (common for PKL PPC models)	
Microplate attachment for one 96 well plate (with base)	
Tube attachment for 20 tubes (with base)	
Round Disk attachment (with base)	
Flask holder attachment (with base)	