

## Digital Tube Rotator

# PKL PPC PR22 Tube Rotator

### Features

- AC synchronous motor for noiseless operation.
- Can accommodate upto 22 Blood collection tubes
- Circular motion gives gentle but effective mixing
- Fixed speed of 30 RPM with an accuracy of +/-2 RPM
- Continuous / infinite Time mode of operation
- Robust construction with strong metal frame
- Plate option 22 x 15 ml or 4 x 50 ml + 12 x 15 ml
- Compact Design - saves bench space
- Protection class of IP 21
- Easy to remove and install plates for replacement and cleaning
- Maximum operational weight capacity of 1kg



### Specifications

| Specifications                                     | Tube Rotator PR22      |
|----------------------------------------------------|------------------------|
| Capacity (supplied with any one chosen disc plate) | 12 x 15 ml + 4 x 50 ml |
| Motor Type                                         | 22 x 15 ml             |
| Motor Type                                         | AC Synchronous         |
| Variable Speed                                     | Fix Speed of 30 RPM    |
| Run Time                                           | Continuous / Infinite  |
| Speed Accuracy                                     | ± 2 RPM                |
| Removable Plate                                    | Yes                    |
| Noise Level                                        | Noiseless operation    |
| Max Loading Capacity)                              | 1 kg                   |
| Dimension (W x D x H)                              | 346 x 185 x 145 mm     |



**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**