

## Overhead Stirrer Series Prime

# **PKL** PPC PPC Prime Overhead Stirrer

### Prime Series Overhead Stirrer

PKL PPC OH 8P (10 liters), OH15P (20 liters), OH 20P (25 liters), OH 40P (40 liters)

#### Features

- Brushless DC motor for quiet operation and maintenance free long operation
- Microprocessor controlled for maintenance of accuracy even with viscous liquid
- Viscosity up to 30,000 mPas & 40 L volume
- Adjustable speed ranging from 50 Up to 2000 RPM
- Long timer setting from 1min to 99 hours 59 mins & continuous mode
- With reversible mixing motion - CW/CCW
- Large digital display for all the parameters at one glance
- Uninterrupted operation through Smart control
- Safety cut for overload, over current & over temperature - safety cut
- Counter reaction to speed change
- Push through shaft for viscous liquid or smaller vessels



**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**

<b>Specifications</b>	<b>Stir OH 08 Prime</b>	<b>Stir OH 15 Prime</b>	<b>Stir OH 20 Prime</b>	<b>Stir OH 40 Prime</b>
Motor Type	Brushless DC	Brushless DC	Brushless DC	Brushless DC
Capacity	10 L	20 L	25 L	40 L
Speed Range	50-2000 RPM	50-1000 RPM	50-2000 RPM	50-1000 RPM
Torque	8 Ncm	15 Ncm	20 Ncm	40 Ncm
Viscosity	4000 mPas	8000 mPas	10000 mPas	30000 mPas
Timer	1 min to 99 hr 59 mins & infinity mode	1 min to 99 hr 59 mins & infinity mode	1 min to 99 hr 59 mins & infinity mode	1 min to 99 hr 59 mins & infinity mode
Direction of Rotation	Reversible - CW/CCW	Reversible - CW/CCW	Reversible - CW/CCW	Reversible - CW/CCW
Protection Class	IP 54	IP 54	IP 54	IP 54
Dimension(Head)	79 x225 x 179 mm	79 x225 x 179 mm	79 x225 x 179 mm	79 x225 x 179 mm
Weight	1.7 kgs	1.7 kgs	1.7 kgs	1.7 kgs

**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**

Via Antonio Amato , 24  
 84131 SALERNO - ITALY  
 Vat.: IT03117920656  
 Ph + 39 089 385027 [amministrazione@paramedical.it](mailto:amministrazione@paramedical.it) [www.paramedical.it](http://www.paramedical.it)

