

## Magnetic Stirrer

# PKL PPC UNO

## Magnetic Stirrer

### Features

- Motorless magnetic stirrer with very slim profile
- No moving parts - Truly maintenance free
- Adjustable speed from 15 to 1500 RPM in steps of 50 RPM
- Timer range 1- 99 mins & infinite mode
- Stirring Capacity of upto 800 ml
- IP 65 compliant design (Water spill resistant)
- Flat surface - Easy to wipe and clean
- Supplied with high strength PTFE Magnetic bar (25 x 8 mm)



### Specifications

Speed Setting

Display

Max. Capacity

Recommended Stir bar length

Pulse Function

Timer

Last run memory

Dimension (W X D X H)

Weight

Protection Class

Standard Accessory

### Magnetic Stirrer UNO

15 to 1500 RPM (at every 50 RPM)

Digital - shows real-time speed

800 ml

25 mm

Yes

1 to 99 mins & infinite mode

Yes

156 x 110 x 13.2 mm

400 gm

IP 65 compliant design (Water spill resistant)

1 high strength stir bar

**PARAMEDICAL srl**

**ISO9001:2015 ISO13485:2016**