


4 Tube Roller

PKL PPC RNR4 Tube Roller

Features

- Gentle rocking and rolling motion for optimum mixing
- Continuous / infinite mode of operation
- Fixed speed of 30 RPM with speed accuracy of +/-2 RPM
- Tilt angle of +/- 3°
- Adjustable rollers spacing to hold different diameter tubes
- Suitable for mixing of blood sample and viscous substances
- Small footprint saving laboratory space
- Works with 80% humidity and 40°C
- Protection class of IP 21
- Robust construction for stable operation



Specifications

Capacity

Motor Type

Variable Speed

Run Time

Tilt angle

Speed Accuracy

Noise Level

Load Capacity

Roller Size (mm)

Protection class

Dimensions (W X D X H)

Weight

Tube Roll RnR 4

9 x 5 ml tubes

AC Synchronous

Fixed Speed of 30 RPM

Continuous / Infinite

+/-30

+/-2 RPM

Noiseless operation

1.5 kg

300 (length) x 30 (roller dia) mm

IP 21

475 x 210 x 80 mm

2 kg

PARAMEDICAL srl

ISO9001:2015 ISO13485:2016

